# Uriline 3 Coli (URILINE 3 COLI)

FR

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GB

Transport medium. Enumeration of microorganisms of the urinary tract, selective growth of enterobacteria and direct identification of *Escherichia coli*.

## SUMMARY AND EXPLANATION

Uriline 3 Coli consists of a closed tube containing a double-sided slide covered with the following culture media:

- Side 1: - CLED agar (green).
- Side 2: - MacConkey agar (reddish brown)  
- *Escherichia coli* detection agar (beige to pale yellow).

This product is intended for use with urine specimens (3) and enables:

- Growth and enumeration of microorganisms of the urinary tract on CLED agar.
- Selective growth and differentiation of enterobacteria on MacConkey agar.
- Direct identification of *Escherichia coli* (2, 4).

## PRINCIPE

### CLED agar:

Also enables lactose-fermenting bacteria to be differentiated from non-fermenting bacteria.

Lactose (+) bacteria produce yellow colonies by acidification of the medium.

Non-lactose fermenting bacteria produce green, blue or colorless colonies.

### MacConkey agar:

Enables detection of lactose fermentation by a change in color of the neutral red indicator.

Lactose-fermenting microorganisms produce red colonies.

Microorganisms that do not ferment lactose produce colonies that are colorless or slightly beige.

The selectivity for Gram (+) bacteria is provided by the bile salts.

### *Escherichia coli* agar:

The presence of a  $\beta$ -glucuronidase substrate enables the identification of *E. coli* colonies: dark brown colonies (3).

The presence of bile salts in the medium inhibits the growth of most Gram (+) bacteria.

## PRESENTATION

REF 56 527 Pack of 10 agar dipslides  
+ 10 blank labels  
+ 1 package insert

## COMPOSITION

Theoretical formula in g/l of purified water.

This medium can be adjusted and/or supplemented according to the performance criteria required:

### CLED agar:

Peptone (bovine or porcine) .....	8
Meat extract (bovine) .....	3
Yeast extract .....	2
Lactose (bovine) .....	10
Cystine .....	0,13
Bromthymol blue .....	0,03
pH 6.9 - 7.4	

### MacConkey agar:

Peptone (bovine or porcine) .....	20
Lactose (bovine) .....	10
Bile salts (bovine) .....	0,8
Neutral red .....	0,075
pH 7.1	

### *E. coli* agar:

Peptone (bovine or porcine) .....	10
Magnesium sulfate .....	0,1
Manganese chloride .....	0,01
Bile salts (bovine) .....	2,4
Ferric citrate .....	0,4
8-Hydroxyquinoline- $\beta$ D-glucuronide .....	q.s.
pH 7.0	

## MATERIAL REQUIRED BUT NOT PROVIDED

- Bacteriology incubator.

## WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to «NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood,*

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*Body Fluids, and Tissue; Approved Guideline – Current Revision*». For further information on handling precautions, refer to «Biosafety in Microbiological and Biomedical Laboratories, CDC/NIH – Latest Edition», or the current regulations in the country of use.

- Culture media should not be used as manufacturing material or components.
- Do not use reagents past the expiration date.
- Do not use contaminated slides, or slides that exude moisture.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient's history, colonial and microscopic morphology and, if necessary, the results of any other tests performed.

## STORAGE CONDITIONS

- **The slides can be stored in their box at 18-25°C until the expiration date.**
- **Avoid variations in temperature and protect from draughts.**

## SPECIMENS

The medium should be directly inoculated using urine (3). Good laboratory practices for collection and transport should be respected.

## INSTRUCTIONS FOR USE

Inoculate the specimen immediately after reception in the laboratory.

1. Unscrew the cap and remove the slide from the tube, avoiding contact with the agar surfaces.
2. Totally immerse the slide in the urine specimen. If the quantity of urine is insufficient, pour the urine over both agar surfaces.
3. Allow excess urine to drain off the slide.
4. Remove the last drops with clean absorbent paper.
5. Insert the slide back into the tube and screw the cap back on tightly.
6. Write the details of the specimen on the label provided and stick it on the tube.
7. Incubate the tube in an upright position at 37°C.  
The user is responsible for choosing the appropriate incubation temperature for the intended use, in accordance with current standards. The cultures are examined after 16-24 hours of incubation.
8. Remove the slide from the tube and perform the reading.

## Storage and transport of inoculated slides:

Prior to incubation, inoculated slides can be stored and transported at 7-25°C. Do not exceed 24 hours following inoculation.


## READING AND INTERPRETATION

After incubation, observe the bacterial growth.

### CLED agar:

#### Enumeration:

Estimate the bacterial concentration by comparing the colony density on the CLED agar surface with the density model chart:

Interpretation of results is as follows:

- Presence of < 10<sup>4</sup> bacteria/ml is normally considered pathologically insignificant.
- Presence of 10<sup>4</sup>-10<sup>5</sup> bacteria/ml is equivocal. Collect another specimen and repeat culture.
- Presence of > 10<sup>5</sup> bacteria/ml indicates probable infection, provided all test procedure instructions have been followed correctly.

#### Aspect of colonies:

- Lactose (+) colonies: yellow.
- Lactose (-) colonies: green, blue or colorless.

### MacConkey agar:

- Record the presence of any colonies.
- Lactose (+) colonies: red.
- Lactose (-) colonies: colorless or slightly beige.

### *E. coli* agar:

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- Record the presence of any characteristic colonies:
- Positive reaction: dark brown colonies (presence of *Escherichia coli*)

- Negative reaction: colorless or beige colonies.

Identification of the microorganism(s) isolated other than *E. coli* must be performed using biochemical or immunological tests.

## QUALITY CONTROL

### Protocol:

The nutrient capacity of the medium can be tested using the following strain:

- *Escherichia coli* ATCC 25922.

### Range of expected results:

Agar	Strain	Results at 33-37°C
<b>CLED</b>	<i>Escherichia coli</i> ATCC 25922	Yellow colonies
<b>MacConkey</b>		Red colonies
<b>E. coli</b>		Dark brown colonies

### Note:

It is the responsibility of the user to perform Quality Control taking into consideration the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature...).

### LIMITATIONS OF THE METHOD

- When the bacterial content is high, the agar surfaces are covered with confluent growth which can easily be overlooked. Consequently, all surfaces which appear sterile when reading is first performed should be examined against reflected light; the absence of reflection indicates confluent growth. This procedure also enables detection of growth of very small colonies.
- MacConkey agar and *E. coli* agar are inhibitory media and must not be used to perform bacterial enumeration.
- Approximately 3% of  $\beta$ -glucuronidase-producing colonies are not *E. coli* (2).
- Approximately 5% of *E. coli* are  $\beta$ -glucuronidase-negative (colorless colonies) (2).
- Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements may not develop.

## PERFORMANCE

### CLED agar:

Detection of bacteriuria at room temperature (1):

Number of samples:	140
Sensitivity:	100%
Specificity:	99%
Reference method:	Petri dish inoculated using a calibrated loop

### E. coli agar:

Performance was evaluated using 2536 strains obtained from urine specimens (2).

### Results obtained:

Microorganisms	Number of strains	Dark brown coloration Number	%
<i>Escherichia coli</i>	1100	1050	95.5
Other enterobacteria	707	10	1.4
Non-fermenting bacilli	284	21	7.4
Gram (+) bacteria	345	0	
Yeasts	100	0	

### WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

### WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

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