**Gélose Trypcase Soja + 5% de sang de mouton (TSS)**

IVD

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**Trypcase Soy Agar + 5% sheep blood (TSS)**

IVD

Isolation of bacteria. Determination of hemolysis

**SUMMARY AND EXPLANATION**

Trypcase soy agar + 5% sheep blood is an isolation medium which favors the growth of all the species currently found in clinical specimens (1, 4).

**PRINCIPLE**

It contains a mixture of peptones particularly suited to the culture of microorganisms. Owing to the presence of sheep blood, hemolysis can be revealed. This is a basic criterion for orienting bacterial identification (2, 3).

**CONTENT OF THE KIT**

Ready-to-use media	
REF 43 001	Pack of 2x10 plates (90 mm)
REF 43 009	Pack of 10x10 plates (90 mm)
<b>TSS *</b>	

\* printed on each plate

**COMPOSITION****Theoretical formula:**

This medium has been adjusted and/or supplemented according to the performance criteria required.

Casein peptone (bovine).....	15 g
Soy peptone .....	5 g
Sodium chloride.....	5 g
Agar .....	15 g
Blood (sheep).....	50 ml
Purified water.....	1 l

pH 7.3

**MATERIAL REQUIRED BUT NOT PROVIDED**

- Controlled atmosphere generators.
- Jars.
- Bacteriology incubator.
- Thermoregulated chambers with a controlled atmosphere.

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI/NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue ; Approved Guideline - Current Revision*".

For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories, CDC/NIH, Latest Edition", or to the regulations currently in use in each country.

- Culture media should not be used as manufacturing material or components.
- Do not use reagents past the expiry date.
- Do not use reagents if the packaging is damaged.
- Do not use contaminated or hemolyzed plates, or plates that exude moisture.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient's history, the source of the specimen, colonial and microscopic morphology and, if necessary, the results of any other tests performed.

**STORAGE CONDITIONS**

- **Store the plates in their box at 2-8°C until the expiry date.**
- If not in the box, plates can be stored for 2 weeks at 2-8°C in the cellophane sachet.

**SPECIMENS**

All types of specimens may be used and should be inoculated directly onto the agar.

Good laboratory practices for collection and transport should be respected and adapted to the type of specimen. This medium can be used to subculture bacterial strains in order to obtain pure cultures.

**INSTRUCTIONS FOR USE**

1. **Allow plates to come to room temperature.**
2. Inoculate the specimen.
3. Put the plate in a suitable atmosphere, if necessary using a controlled atmosphere generator.
4. Incubate with the cover bottom side at 37°C. The user is responsible for choosing the appropriate incubation temperature for the intended use, in accordance with current standards. Incubation time varies according to the type of specimen and the microorganisms being tested for. The cultures are generally examined after 24-48 hours of incubation. In certain cases, it may be necessary to prolong incubation.

### READING AND INTERPRETATION

- After incubation, observe the bacterial growth.
- Record the presence of any characteristic hemolyses:
  - $\alpha$  hemolysis: greenish discoloration around the colony.
  - $\beta$  hemolysis: clear zone around or under the colony.
- Identification of the microorganism(s) isolated must be followed by biochemical and/or immunological tests.

### QUALITY CONTROL

#### Protocol:

The nutrient capacity of the medium can be tested using the following strains (incubation in a CO<sub>2</sub>-enriched atmosphere):

- *Streptococcus pyogenes* ATCC® 19615
- *Streptococcus pneumoniae* ATCC® 6305

#### Range of expected results:

Strain	Results at 33-37°C	
<i>Streptococcus pyogenes</i> ATCC® 19615	Growth after 24 hours	$\beta$ hemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305		$\alpha$ hemolysis

#### Note:

It is the responsibility of the user to perform Quality Control taking into consideration the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature, etc.).

### LIMITATIONS OF THE METHOD

- Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements (substrate, temperature, incubation conditions, etc.) may not develop.
- The type of hemolysis depends on the species in question and the specific behavior of each strain.
- Depending on the specimens analyzed and the microorganisms being tested for, it is recommended to use Trypcase soy agar + 5% sheep blood in conjunction with additional media (selective media, Chocolate agars, etc.).

### PERFORMANCE

Performance was evaluated at 37°C using 34 bacterial strains (*Streptococcus pyogenes*, *Streptococcus agalactiae*, group C and D streptococci, *Streptococcus pneumoniae*, staphylococci and other Gram (+) cocci, Gram (-) bacilli, *Listeria*, anaerobic bacteria) and 1 yeast (*Candida*).

### Nutrient capacity:

All the strains grew within 24-48 hours.

### Hemolysis:

The following types of hemolyses are observed with strains tested after 24 hours of incubation:

- $\beta$  hemolysis for the 5 strains of *Streptococcus pyogenes* and the 4 strains of *Listeria* (*L. monocytogenes*, *L. seeligeri*, *L. ivanovii*).  
Note:  $\beta$  hemolysis is not characteristic of these species: certain staphylococci, *Streptococcus agalactiae* and group C streptococci have also produced  $\beta$  hemolysis.
- $\alpha$  hemolysis for the 5 strains of *Streptococcus pneumoniae*.

### WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

### LITERATURE REFERENCES

- 1 ATLAS R.M., *Handbook of microbiological media* – CRC Press, 1993.
- 2 FACKLAM R.R., PADULA J.F., WORTHAM E.C. and al. – Presumptive identification of group A, B, and D streptococci on agar plate media. – *J. Clin. Microbiol.*, 1979, vol. 9, n°6, p. 665-672.
- 3 FLANDROIS J.P., CHOMARAT M. – *Bactériologie médicale pratique* - MEDSI / Mac GRAW-HILL, 1989 – ISBN 2-86439-161-9.
- 4 MURRAY P.R., BARON E.J., PFALLER M.A. et al. – *Manual of clinical microbiology* – 6th ed. – Washington D.C. : ASM Press, 1995 – ISBN 1-55581-086-1.

**INDEX OF SYMBOLS**

Symbol	Meaning
 or REF	GB : Catalogue number US : Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests

**WARRANTY**

*bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.*

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

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bioMérieux and the blue logo are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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