Gélose Trypcase Soja (TSA-T)

IVD

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Trypcase Soy Agar (TSA-T)

IVD

Culture of nonfastidious microorganisms

SUMMARY AND EXPLANATION

Trypcase Soy agar is a medium for facilitating the growth of microorganisms which do not have specific growth requirements (1).

This medium corresponds to the B agar medium of the European Pharmacopoeia (4). It also complies with the USP (5).

PRINCIPLE

It contains a mixture of peptones which enable the growth of most microorganisms.

CONTENT OF THE KIT

	Ready-to-use medium
REF 42 101	20 tubes (slanted)

COMPOSITION

Theoretical formula.

This medium can be adjusted and/or supplemented according to the performance criteria required:

Casein peptone (bovine).....	15 g
Soy peptone	5 g
Sodium chloride.....	5 g
Agar.....	15 g
Purified water.....	1 l

pH 7.3

MATERIAL REQUIRED BUT NOT PROVIDED

- Controlled atmosphere generators.
- Jars.
- Bacteriology incubator.

Or

- Thermoregulated chambers with a controlled atmosphere.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI/NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline - Current Revision*". For additional information on handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories – CDC/NIH – Latest edition", or the current regulations in the country of use.

- Culture media should not be used as manufacturing material or components.
- Do not use reagents past the expiry date.
- Do not use tubes which show signs of contamination.
- Before use, check that the cap is intact.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient's history, the source of the specimen, colonial and microscopic morphology and, if necessary, the results of any other tests performed.

STORAGE CONDITIONS

- **Store the tubes in their box at 2-25°C until the expiry date.**

SPECIMENS

This medium may be used to subculture microbial strains.

INSTRUCTIONS FOR USE

1. Inoculate the tube with the strain to be cultured.
2. Put the tube in a suitable atmosphere, if necessary using controlled atmosphere generators.
3. Incubate at 37°C. The user is responsible for choosing the appropriate temperature for the intended use, in accordance with current standards.
Incubation time varies according to the type of microorganism. The cultures are generally examined after 24-48 hours of incubation. In certain cases, however, incubation may be prolonged.

READING AND INTERPRETATION

After incubation, observe the growth.

QUALITY CONTROL

Protocol:

The nutrient capacity of the medium can be tested using the following strain:

Staphylococcus aureus ATCC® 6538

Range of expected results:

At 33-37°C, the tested strain should grow after 24 hours of incubation.

Note:

It is the responsibility of the user to perform Quality Control taking into consideration the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature, etc.).

LIMITATIONS OF THE METHOD

Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements may not develop.

PERFORMANCE

Performance was evaluated at 37°C using 16 bacterial strains: 7 Gram (-) strains (enterobacteria and other Gram (-) organisms), 9 Gram (+) strains, and 1 yeast (*Candida*).

Nutrient capacity:

All the tested strains grew within 24 hours.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

1. FINEGOLD S.M., MARTIN W.J., SCOTT E.G. – *Bailey and Scott's Diagnostic microbiology* – MOSBY – 1978.
2. Mac FADDIN J.F. - 1985 - *Media for isolation-cultivation-identification-maintenance of medical bacteria*, vol. 1, Williams and Wilkins, Baltimore.
3. MURRAY P.R., BARON E.J., PFALLER M.A. et al. - 1995 – *Manual of clinical microbiology, 6th ed.* - American Society for Microbiology – ISBN 1-55581-086-1.
4. European Pharmacopoeia IV (2002) – ISBN 92-871-4841-4.
5. USP XXV (2002).

INDEX OF SYMBOLS

Symbol	Meaning
	GB : Catalogue number US : Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

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bioMérieux and the blue logo are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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