




Wir sind für Sie da:  
E-Mail: [office@rieger-iv.at](mailto:office@rieger-iv.at)  
Telefon: +43(1) 728 00 52  
Fax: +43(1) 728 69 16  
Rustenschacher Allee 10  
1020 Wien, Austria



# DISTRIWEL 440<sup>®</sup>

## Media Pourer Stacker





# DISTRIWEL 440®

Increase your productivity  
quality  
traceability

## PERFORMANCES



440 plates  
poured  
in 30 mins

# DISTRIWEL 440®

### Walkaway device!

Up to 440 dishes  
filled in 30 mins



0°



Ensure a flat poured media  
(bubble level integrated)

7" touch screen LCD  
Intuitive display

**Traceability**  
Print up to 25 characters  
(2 lines) on the Petri dish  
(optional laser printer)



The Built-in Peltier cooling  
effect reduces condensation  
and solidification time

### AEROBIOCONTAMINATION PREVENTION

Contamination prevention  
(UV lamp) while pouring  
+ Plexiglas protection cover

In-line blood  
addition option


Multicolour  
visual alarms

Your culture media is accurately  
poured by a Masterflex™  
peristaltic pump



## SPECIFICATIONS

# DISTRIWEL 440®



## DISTRIWEL FEATURES

Dishes/hour for 90mm	Up to 750
Capacity of the carousel	440 dishes
Fast plate cooling	Built-in Peltier effect
Contamination prevention	Built-in UV lamp
Pour plates function	✓
Pouring accuracy	±1%
Pouring range	1 – 99 ml
Stack plates function	✓
Peristaltic pump	Integrated Masterflex pump
Stored programs	40
Control panel	LCD Touch screen
Country of origin	100% France
Traceability (optional)	Laser printer (LINX CJ400) → Printed Petri dishes
Footswitch (for manual dispensing)	✓ (optional)

## DIMENSIONS & ELECTRICAL SUPPLY

W x H x D	808 x 940 x 689 mm
Electrical supply	100-250v 50-60 Hz 10 A
Power supply	150 W

PACK YOUR  
READY  
TO USE DISHES !



ORDER IT TOGETHER WITH DISTRIWEL 440 : Ref. MINIMA



## ALLIANCE BIO EXPERTISE



Empowered from our people's skills and experiences, we, ALLIANCE BIO EXPERTISE, together with our customers and partners decided to manufacture Microbiology automated equipment around essential values:  
Customer Focus, Service, Loyalty, Innovations

**Ask your local partner  
for a personal demonstration!**

ZI de courbouton , Le Tremplin – 35480 GUIPRY - FRANCE  
E-mail : [contact@abioexpertise.com](mailto:contact@abioexpertise.com)  
Phone : +33(0)2.40.51.79.53  
[alliance-bio-expertise.com](http://alliance-bio-expertise.com)