

EasyCount 2

EC2™

*The color colony counter that saves time,
bench space and money!*



- Accurate and rapid Counting (± 1 sec)
- Cost-effective
- Full traceability (results, parameters, image)
- Instant report generation (spreadsheet and LIMS compatible)
- Ideal for pour plates and spiral plates
- Compact footprint saves bench space


BIOMÉRIEUX
INDUSTRY

PIONEERING DIAGNOSTICS

EasyCount 2 standardizes your colony counting

• **Meets your needs for:**

- Fully automatic or semi-automatic counting
- Versatility - Ideal for pour plates and spiral plates

• **Improve quality and productivity:**

- Provides automatic accurate counting quickly (± 1 sec per plate)
- Ability to perform a whole plate or sector count

• **Experience total ease of use:**

- Plug & play installation through USB connectivity
- Live, full color on-screen images of the counted colonies
- Plate image zoom available
- Interactive software
- Selectable top or base illumination (upper and lower LED arrays with darkfield and brightfield)



Traceability

Automatically save results, plate identification and plate image

Optional bar code reader for instant reading of plate information

Result reports can be printed and are available for further data processing (spreadsheet or LIMS)



Ordering Information

AESEM1011

Automatic Colony Counter EC2

Technical specifications

For 55 or 90 mm diameter plates

Power: 110/240 V

Weight: 7 lbs

Dimensions: 10" x 11" x 10"

Minimum PC configuration: 500 MHz Pentium III, RAM MB128, CD-ROM drive, Windows 2000, XP

2 free USB ports, 1024 x 768, 24 bit true color graphics display

bioMérieux, Inc.

595 Anglum Road
Hazelwood, MO 63042 • U.S.A.
Tel: (800) 634 7656
Fax: (800) 657 3053

www.biomerieux-usa.com

bioMérieux Canada, Inc.

7815 Henri-Bourassa Boulevard West
St. Laurent (Québec) Canada H4S 1P7
Tel: (1) 514 336 7321
Fax: (1) 514 807 0015

